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ENGLISH

INTENDED USE

RNA Medical® Brand QC 463 EQUIL Plus with Glucose and Lactate Control is a quality control material that, once tonometered, is used for monitoring the performance of blood gas, electrolyte, metabolite, and CO-Oximeter instrumentation for the analytes and analyzers listed on the Expected Values Chart.

PRODUCT DESCRIPTION

QC 463 is provided in three (3) levels for monitoring analyzer performance at different points within the clinical range. It is packaged in sealed glass ampuls, each containing 2.4 mL of solution. Ampuls are packaged thirty (30) per box.

Active Ingredients:

QC 463 is a buffered bovine hemoglobin solution containing electrolytes (Na⁺, K⁺, Cl⁻, Ca⁺⁺, Mg⁺⁺), glucose, and lactate. Prior to use, QC 463 is tonometered with precision gas mixtures of CO₂, O₂, and N₂ to produce distinct control levels for pH, *p*CO₂, *p*O₂, and O₂Hb. This control contains no preservatives or human-based materials. It is considered good laboratory practice to follow the recommended "Universal Precautions" when handling any blood product.

STORAGE

The expiration date stated on the QC 463 packaging is for product stored refrigerated (2-8 °C). The product may remain at room temperature (up to 25 °C) for seven (7) days, provided the labeled expiration date is not exceeded. Avoid exposure to freezing and prolonged exposure to temperatures greater than 8 °C.

MATERIALS REQUIRED BUT NOT PROVIDED

QC 463 is for use in tonometry. The RNA Medical EQUILibrator® Tonometer is designed for use with EQUIL controls. It requires compressed gas cylinders with precision (+/- 0.1%) gas mixtures, suitable regulators, and several other supplies for proper operation. Refer to the EQUILibrator Operator's Manual for a detailed account of required accessories.

DIRECTIONS FOR USE

The procedure below is for use of QC 463 with the EQUILibrator Tonometer. Detailed instructions for tonometry of this and all EQUIL controls are contained in the EQUILibrator Operator's Manual.

Prior to tonometry, the control should be brought to room temperature (18-25 °C). Allow at least two (2) hours for the ampuls to equilibrate to this temperature.

Using gauze, tissue, gloves, or an appropriate ampul opener to protect fingers from cuts, open the ampul by snapping off the tip at the score. Aspirate the contents into a Pre-Treated Syringe (RNA Catalog # T 310) and tonometer the sample as directed. After tonometry, introduce the prepared control directly from the syringe to the analyzer, following the instrument manufacturer's instructions for sampling a control material.

All parameters (including *p*O₂) of tonometered QC 463 will remain stable for up to one (1) hour in a capped syringe. Less volume in the syringe will result in less stability.

EXPECTED VALUES

The blood gas values for this control are calculated as follows:

*p*O₂ (or *p*CO₂) in mmHg = [ambient barometric pressure - 47 (the partial pressure of water vapor at 37 °C)] X % O₂ (or % CO₂) in the gas mixture

The values for all other control analytes on the enclosed Expected Values Chart are based on multiple determinations performed on randomly selected samples from each lot. The listing for each instrument represents the expected range and mean value of this range.

The Expected Values are provided as a guide in evaluating analyzer performance. Since instrument design and operating conditions may vary, each laboratory should establish its own expected values and control limits. The mean value established should fall within the Expected Value range shown on the chart.

RNA Medical provides monthly statistical reports for tracking and review of analyzer performance as well as lot number specific peer group data. Please contact RNA Medical for information about this service.

LIMITATIONS

- QC 463 is sensitive to many instrument related factors that would affect analytical results. Although it is a bovine blood-based material, it does not contain red cells. Therefore, it may not detect certain malfunctions that would affect the testing of human blood.
- This product is intended for use as a quality control material and can assist in evaluating the performance of laboratory instruments. It is not for use as a calibration standard and its use should not replace other aspects of a complete quality control program.
- To avoid foaming, QC 463 should be equilibrated in RNA Medical Pre-Treated Syringes (Catalog # T 310).
- Accurate values will only be obtained when QC 463 has been tonometered for the recommended time and at the correct temperature without dilution or concentration of the control.

RNA Medical, EQUIL, and the EQUILibrator are registered trademarks of Bionostics, Inc. The products described herein are covered by one or more of the following U.S. Patents and their foreign counterparts: 5,558,985; 5,320,965; 5,304,491; 5,945,529; 5,913,866; 4,945,062.

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