QC 253 Full Range CO-Oximeter Control

ESPAÑOL

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ENGLISH

INTENDED USE

RNA Medical® Brand QC 253 Full Range CO-Oximeter Control is an assayed quality control material used for monitoring the measurements of total hemoglobin and hemoglobin fractions on CO-Oximeter instrumentation for the analyzers listed on the Expected Values Chart.

PRODUCT DESCRIPTION

QC 253 is provided in three (3) levels of total hemoglobin, oxyhemoglobin, and carboxyhemoglobin for monitoring analyzer performance at different points within the clinical range; it also contains methemoglobin. QC 253 is packaged in sealed glass ampuls, each containing 1.2 mL of solution. Ampuls are packaged thirty (30) per box.

Active Ingredients:

QC 253 is a purified bovine hemoglobin solution which has been treated with compressed air and precise concentrations of carbon monoxide. This control contains no preservatives and no human-based materials. It is considered good laboratory practice to follow the recommended "Universal Precautions" when handling any blood product.

STORAGE

The expiration date stated on the QC 253 packaging is for product stored refrigerated (2- 8 °C). Avoid exposure to freezing and temperatures greater than 8 °C.

DIRECTIONS FOR USE

QC 253 should be analyzed immediately after removal from refrigeration. Before use, gently invert the ampul to mix the solution. Tap the ampul to restore the liquid to the bottom. Use gauze, tissue, gloves, or an appropriate ampul opener to protect fingers from cuts and open the ampul by snapping off the tip at the score. Immediately introduce the liquid from the ampul to the analyzer, following the instrument manufacturer's instructions for sampling a control material. Use direct aspiration, syringe transfer, or capillary mode techniques.

EXPECTED VALUES

The values for each control analyte on the enclosed Expected Values Chart are based on multiple determinations performed on randomly selected samples from each lot. The listing for each instrument represents the expected range and mean value of this range.

The Expected Values are provided as a guide in evaluating analyzer performance. Since instrument design and operating conditions may vary, each laboratory should establish its own expected values and control limits. The mean value established should fall within the Expected Value range shown on the chart.

STATISTICAL SUPPORT

RNA Medical PeerQC®, available at www.RNAMedical.com, provides monthly statistical reports for tracking and review of analyzer performance as well as lot number specific peer group data. This service is available at no charge to RNA Medical customers. Please contact RNA Medical or visit our website for information about utilizing PeerQC for this product.

LIMITATIONS

- QC 253 is sensitive to many instrument related factors that would affect analytical results. Although it is a bovine blood-based material, it does not contain red cells. Therefore, it may not detect certain malfunctions that would affect the testing of human blood.
- This product is intended for use as a quality control material and can assist in evaluating the performance of laboratory instruments. It is not for use as a calibration standard and its use should not replace other aspects of a complete quality control program.
- Extended exposure to temperatures greater than 8 °C will affect product performance. If QC 253 has turned brown in color, this change indicates deterioration and the formation of methemoglobin. In such a case, the control is unsuitable for use and should be discarded.

RNA Medical is a registered trademark and PeerQC is a registered service mark of Bionostics, Inc., Devens, MA, USA. The products and services described are covered by U.S. Patent: 7,027,917

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Siemens	tHb g/dL	O ₂ Hb %	COHb %	MetHb %	HHb %
400 Series	14.3	13.3 - 15.3	83.4	79.4 - 87.4	16.6
500 Series	DNA ¹	DNA ¹	DNA ¹	DNA ¹	DNA ¹
1200 Series	14.2	13.2 - 15.2	83.4	79.4 - 87.4	16.6

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Roche	tHb g/dL	COHb %	MetHb %	O ₂ Hb %	RHb (HHb) %
OMNI Series	16.2	15.0 - 17.4	43.4	39.4 - 47.4	1.1
				-0.9 - 3.1	55.0
				51.0 - 59.0	DNA ¹

Siemens	tHb g/dL	O ₂ Hb %	COHb %	MetHb %	HHb %
400 Series	17.4	16.2 - 18.6	56.5	52.5 - 60.5	44.4
500 Series	DNA ¹	DNA ¹	DNA ¹	DNA ¹	DNA ¹
1200 Series	17.5	16.3 - 18.7	56.5	52.5 - 60.5	44.4

Footnotes / Anmerkungen / Notas al pie / Notes explicatives / Note / Fotnoter:

1. Data not available at the time of printing. / Daten bei Drucklegung nicht vorhanden. / Datos no disponibles en el momento de la impresión. / Aucune donnée disponible lors de cette publication. / Dati non disponibili al momento della stampa. / Uppgifter fanns ej tillgängliga vid tidpunkten för tryckning.