



May 8, 2017

Dear RNA Medical Customer,

Please find the AADEE Expected Values for RNA Medical QC123 products below.

	AADEE Rumi BG/AADEE ZEN								
	pH			pCO2 mmHg			pO2 mmHg		
	Media Mean	Rango Range		Media Mean	Media Mean		Rango Range	Media Mean	
QC123-1 Lot #51815 Exp: 2019-06	7,26	7,23	7,29	60	53	67	75	65	85
QC123-2 Lot #51915 Exp: 2019-06	7,47	7,43	7,50	38	31	45	114	104	124
QC123-3 Lot #52016 Exp: 2019-06	7,67	7,63	7,70	23	16	30	144	134	154

Please do not hesitate to contact RNA Medical with any questions:

RNA Medical, Division of Bionostics, Inc.
 7 Jackson Road
 Devens, MA 01434
 Phone: 978-772-9070 / 800-533-6162
 Fax: 978-598-7918
 E-Mail: customerservice@rnamedical.com