

May 8, 2017

Dear RNA Medical Customer,

Please find the AADEE Expected Values for RNA Medical QC623 products below.

AADEE μ ISE												
	Na+ mmol/L			K+ mmol/L			Cl- mmol/L			Ca++ mmol/L		
	Media Mean	Rango Range		Media Mean	Rango Range		Media Mean	Rango Range		Media Mean	Rango Range	
QC623-1 Lot #40222 Exp: 2018-07	116	111	121	2,2	1,7	2,7	80	75	85	1,52	1,37	1,67
QC623-2 Lot #40322 Exp: 2018-07	134	129	139	4,3	3,8	4,8	93	88	98	1,15	1,05	1,25
QC623-3 Lot #40422 Exp: 2018-08	156	151	161	6	5,5	6,5	123	118	128	0,52	0,42	0,62
QC623-1 Lot #50224 Exp: 2019-06	115	110	120	2,2	1,7	2,7	78	73	83	1,52	1,42	1,62
QC623-2 Lot #50323 Exp: 2019-05	136	129	141	4,2	3,7	4,7	93	88	98	1,15	1,05	1,25
QC623-3 Lot #50423 Exp: 2019-05	160	155	165	5,9	5,4	6,4	123	118	128	0,51	0,41	0,61

AADEE Mini ISE												
	Na+ mmol/L			K+ mmol/L			Cl- mmol/L			Ca++ mmol/L		
	Media Mean	Rango Range		Media Mean	Rango Range		Media Mean	Rango Range		Media Mean	Rango Range	
QC623-1 Lot #50224 Exp: 2019-06	119	114	124	2	1,5	2,5	80	75	85	1,43	1,28	1,58
QC623-2 Lot #50323 Exp: 2019-05	133	128	138	4,2	3,7	4,7	93	88	98	1,18	1,08	1,28
QC623-3 Lot #50423 Exp: 2019-05	151	146	156	6	5,5	6,5	116	111	121	0,6	0,50	0,70

AADEE Rumi BG/ AADEE ZEN									
	pH			pCO2 mmHg			pO2 mmHg		
	Media Mean	Rango Range		Media Mean	Rango Range		Media Mean	Rango Range	
QC623-1 Lot #50224 Exp: 2019-06	7,15	7,12	7,18	65	58	73	79	69	89
QC623-2 Lot #50323 Exp: 2019-05	7,40	7,36	7,43	45	40	50	112	102	122
QC623-3 Lot #50423 Exp: 2019-05	7,60	7,57	7,64	24	19	29	142	132	152

Please do not hesitate to contact RNA Medical with any questions:

RNA Medical, Division of Bionostics, Inc.
 7 Jackson Road
 Devens, MA 01434
 Phone: 978-772-9070 / 800-533-6162
 Fax: 978-598-7918
 E-Mail: customerservice@rnamedical.com