

August 1, 2014

Dear RNA Medical Customer,

AADEE Rumi BG ranges have been established for our QC123 lot# 31814/BR31216 for L1, Lot# 31914/BR31217 for L2 and Lot# 32015/BR31218 for L3. Please see the chart below.

Analyte	LEVEL 1 MEAN	LEVEL 1 RANGE	LEVEL 2 MEAN	LEVEL 2 RANGE	LEVEL 3 MEAN	LEVEL 3 RANGE
pH	7.23	7.20-7.26	7.44	7.41-7.47	7.64	7.61-7.67
PCO ₂ mmHg	57	50-64	38	33-43	23	18-28
PO ₂ mmHg	678	68-88	110	100-120	143	133-153

RNA Medical has not reprinted the insert for this lot, so please use this letter as notification for this update. Please keep this announcement with your records as long as you have these lots in stock.

Please do not hesitate to contact RNA Medical with any questions regarding this update at 1-800-533-6162 or CustomerService@RNAMedical.com.

Sincerely,

Jason Mack

RNA Product Manager

1 de Agosto de 2014,

Estimado RNA Cliente Médico,

Se han establecido rangos AADEE Rumi BG para nuestro QC123 lot# 31814/BR31216 para L1, Lot# 31914/BR31217 para L2 y Lot# 32015/BR31218 para L3. Por favor, consulte la siguiente tabla.

Analyte	LEVEL 1 MEDIA	LEVEL 1 RANGO	LEVEL 2 MEDIA	LEVEL 2 RANGO	LEVEL 3 MEDIA	LEVEL 3 RANGO
pH	7.23	7.20-7.26	7.44	7.41-7.47	7.64	7.61-7.67
PCO ₂ mmHg	57	50-64	38	33-43	23	18-28
PO ₂ mmHg	678	68-88	110	100-120	143	133-153

RNA Medical no ha reproducido el insert para este lote, así que favor de utilizar esta carta como notificación para esta actualización. Por favor mantenga este anuncio con sus registros mientras tienes estos lotes en stock. Por favor, no dude en contactar con RNA Medical con cualquier pregunta relacionada con esta actualización al 1-800-533-6162 o CustomerService@RNAMedical.com.

Sinceramente,

Jason Mack

RNA Product Manager