

**CC 527 Blood Gas/Electrolyte/CO-Oximeter Control (Accriva AVOXimeter and LifeHealth IRMA TRUpoint)**

Level 1

43329

2027-08-31

Expected Values Chart

Manufacturer/ CO-Oximeter	tHb g/dL			O <sub>2</sub> Hb %			COHb %			MetHb %		
	mean	range		mean	range		mean	range		mean	range	
<b>Accriva</b>												
AVOXimeter 4000	6.5	5.5 - 7.4		82.5	78.4 - 86.6		16.1	12.1 - 20.1		2.7	-0.2 - 5.6	
AVOXimeter 1000E	6.8	5.7 - 7.8		83.7	80.2 - 87.2							

Level 2

43429

2027-08-31

Expected Values Chart

Manufacturer/ CO-Oximeter	tHb g/dL			O <sub>2</sub> Hb %			COHb %			MetHb %		
	mean	range		mean	range		mean	range		mean	range	
<b>Accriva</b>												
AVOXimeter 4000	11.2	10.0 - 12.3		88.6	84.9 - 92.2		8.0	4.3 - 11.6		3.5	0.8 - 6.2	
AVOXimeter 1000E	11.5	10.3 - 12.7		88.9	85.3 - 92.4							

Level 3

43529

2027-08-31

Expected Values Chart

Manufacturer/ CO-Oximeter	tHb g/dL			O <sub>2</sub> Hb %			COHb %			MetHb %		
	mean	range		mean	range		mean	range		mean	range	
<b>Accriva</b>												
AVOXimeter 4000	17.7	16.0 - 19.3		94.2	90.9 - 97.5		1.4	-1.5 - 4.3		4.3	1.8 - 6.8	
AVOXimeter 1000E	17.5	15.8 - 19.1		94.2	90.9 - 97.5							

Level 1,2,3

Manufacturer/ Analyzer	pH		pCO <sub>2</sub> mmHg		pO <sub>2</sub> mmHg		Ca <sup>++</sup> mmol/L		Na <sup>+</sup> mmol/L		K <sup>+</sup> mmol/L		Cl <sup>-</sup> mmol/L		H <sup>+</sup> nmol/L		pCO <sub>2</sub> kPa		pO <sub>2</sub> kPa	
	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range
<b>LifeHealth</b>																				
IRMA TRUpoint <sup>1</sup>																				

1. Contact LifeHealth Technical Support or the LifeHealth website (www.lifehealthmed.com) for assigned values. / Die zugeordneten Werte erfahren Sie vom LifeHealth Kundendienst oder von der LifeHealth Website (www.lifehealthmed.com). / Para obtener información sobre los valores asignados, póngase en contacto con el servicio de asistencia técnica de LifeHealth o visite la web de LifeHealth en Internet (www.lifehealthmed.com). / Contacter l'assistance technique de LifeHealth ou son site web (www.lifehealthmed.com) pour les valeurs attribuées. / Per i valori assegnati, rivolgersi al servizio di assistenza tecnica di LifeHealth o visitare il sito www.lifehealthmed.com. / Kontakta LifeHealth tekniska support eller besök LifeHealth webbplats (www.lifehealthmed.com) för information om angivna värden.