

CC 527 Blood Gas/Electrolyte/CO-Oximeter Control (Accriva AVOXimeter and LifeHealth IRMA TRUpoint)

Level 1

13327

2024-03-31

Expected Values Chart

Manufacturer/ CO-Oximeter	tHb g/dL		O2Hb %		COHb %		MetHb %								
	mean	range	mean	range	mean	range	mean	range							
Accriva															
AVOXimeter 4000	6.9	5.9 - 7.8	82.3	78.2 - 86.4	15.6	11.6 - 19.6	3.1	0.2 - 6.0							
AVOXimeter 1000E	7.1	6.0 - 8.1	83.6	80.1 - 87.1											

Level 2

13427

2024-03-31

Expected Values Chart

Manufacturer/ CO-Oximeter	tHb g/dL		O2Hb %		COHb %		MetHb %								
	mean	range	mean	range	mean	range	mean	range							
Accriva															
AVOXimeter 4000	11.5	10.3 - 12.6	88.2	84.5 - 91.8	7.8	4.1 - 11.4	3.4	0.7 - 6.1							
AVOXimeter 1000E	11.8	10.6 - 13.0	88.7	85.1 - 92.2											

Level 3

13527

2024-03-31

Expected Values Chart

Manufacturer/ CO-Oximeter	tHb g/dL		O2Hb %		COHb %		MetHb %								
	mean	range	mean	range	mean	range	mean	range							
Accriva															
AVOXimeter 4000	18.0	16.3 - 19.6	94.3	91.0 - 97.6	0.7	-2.2 - 3.6	4.7	2.2 - 7.2							
AVOXimeter 1000E	18.1	16.4 - 19.7	94.1	90.8 - 97.4											

Level 1,2,3

Manufacturer/ Analyzer	pH		pCO2 mmHg		pO2 mmHg		Ca++ mmol/L		Na+ mmol/L		K+ mmol/L		Cl- mmol/L		H + nmol/L		pCO2 kPa		pO2 kPa		
	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	mean	range	
LifeHealth																					
IRMA TRUpoint ¹																					

1. Contact LifeHealth Technical Support or the LifeHealth website (www.lifehealthmed.com) for assigned values. / Die zugeordneten Werte erfahren Sie vom LifeHealth Kundendienst oder von der LifeHealth Website (www.lifehealthmed.com). / Para obtener información sobre los valores asignados, póngase en contacto con el servicio de asistencia técnica de LifeHealth o visite la web de LifeHealth en Internet (www.lifehealthmed.com). / Contacter l'assistance technique de LifeHealth ou son site web (www.lifehealthmed.com) pour les valeurs attribuées. / Per i valori assegnati, rivolgersi al servizio di assistenza tecnica di LifeHealth o visitare il sito www.lifehealthmed.com. / Kontakta LifeHealth tekniska support eller besök LifeHealth webbplats (www.lifehealthmed.com) för information om angivna värden.