



February 2, 2026

Dear Valued Customer,

QC253 multi-level set lot 548074, Expiration date 2026-10-31 and QC2531 level 1 set lot 54880, Expiration date 2026-10-31 has incorrect target values and ranges for the following Hemoglobin targets and ranges.

Avoximeter 4000, Level 1 O2Hb and COHb, Lot 54880, Expiration date 2026-10-31

|         | Revised Values O2Hb and COHb % |      |       |  | Original Values O2Hb and COHb % |      |       |
|---------|--------------------------------|------|-------|--|---------------------------------|------|-------|
| Analyte | Mean                           | Low  | High  |  | Mean                            | Low  | High  |
| O2Hb    | 96.1                           | 91.6 | 100.6 |  | 97.4                            | 92.9 | 101.9 |
| COHb    | 5.5                            | 1.5  | 9.5   |  | 4.3                             | 0.3  | 8.3   |

The insert that comes with the kit will not contain this information. This letter will serve as notification for the value assignment of these analytes for this lot. Please keep this announcement with your records for as long as you have this lot in stock.

If you have any questions regarding this notification, please feel free to call or email us using the contact information provided.

RNA Medical, Division of Bionostics, Inc.  
7 Jackson Road  
Devens, MA 01434  
Phone: 978-772-9070 / 800-533-6162  
E-Mail: [customerservice@rnamedical.com](mailto:customerservice@rnamedical.com)

BIOMKT-243 Rev. A