



December 09, 2025

Dear Valued Customer,

CVC223 set lot 522345, Expiration date 2026-11-30 has incorrect target values and ranges for level 4 MetHb, lot 55075, Expiration date 2027-01-31. Please find the updated ranges for MetHb in the table below.

Instrument	Revised Set lot 522345 Level 4 lot 55075 MetHb %			Original set lot 522345 Level 4 lot 55075 MetHb %		
	Mean	Low	High	Mean	Low	High
682	0.7	-1.3	2.7	-4.4	-6.4	-2.4
Synthesis Series	0.6	-1.4	2.6	-4.5	-6.5	-2.5
GEM OPL	0.0	-3.0	3.0	0.0	-3.0	3.0
GEM Premier 4000	1.4	-0.6	3.4	-3.7	-5.7	-1.7
GEM Premier 5000	1.4	-0.6	3.4	-3.7	-5.7	-1.7
ABL 700 Series	3.0	1.0	5.0	-2.1	-4.1	-0.1
ABL 800 Series	3.3	1.3	5.3	-1.8	-3.8	0.2
ABL 90 Series	2.2	0.2	4.2	-2.9	-4.9	-0.9
Cobas b 221	1.0	-1.0	3.0	-4.1	-6.1	-2.1
OMNI Series	1.6	-0.4	3.6	-3.5	-5.5	-1.5
400 Series	0.8	-1.2	2.8	-4.3	-6.3	-2.3
500 Series	1.1	-0.9	3.1	-4.0	-6.0	-2.0
1200 Series	1.2	-0.8	3.2	-3.9	-5.9	-1.9

RNA Medical, Division of Bionostics, Inc.
 7 Jackson Road
 Devens, MA 01434
 Phone: 978-772-9070 / 800-533-6162
 E-Mail: customerservice@rnamedical.com

BIOMKT-241 Rev. A