



March 4, 2021

Dear Valued Customer,

AADEE has released the following information for their instrument values for RNA Medical QC623 Kit lot 902776. Please find the expected values and ranges for the QC623 lots below.

AADEE μISE									
	QC623 Lot#	Na ⁺		K ⁺		Cl ⁻		Ca ⁺⁺	
		Mean	Range	Mean	Range	Mean	Range	Mean	Range
Level 1	90227	118	112-124	2.3	1.8-2.8	79	74-84	1.46	1.33-1.63
Level 2	90327	135	130-140	4.3	3.8-4.8	94	86-102	1.12	0.97-1.22
Level 3	90426	159	154-164	6.2	5.7-6.7	122	117-127	0.53	0.42-0.63

AADEE RUMI BG							
	QC623 Lot#	pH		pCO ₂ (mmHg)		pO ₂ (mmHg)	
		Mean	Range	Mean	Range	Mean	Range
Level 1	90227	7.155	7.12-7.19	67	61-73	77	70-87
Level 2	90327	7.40	7.36-7.45	45	39-49	112	103-122
Level 3	90426	7.61	7.57-7.65	26	23-29	147	132-157

AADEE μGases							
	QC623 Lot#	pH		pCO ₂ (mmHg)		pO ₂ (mmHg)	
		Mean	Range	Mean	Range	Mean	Range
Level 1	90227	7.157	7.12-7.19	71	65-76	75	66-86
Level 2	90327	7.41	7.37-7.46	46	41-51	106	98-120
Level 3	90426	7.62	7.57-7.65	21	17-27	145	134-155

No other instruments are affected by this update for QC623 lot 902776.

If you have any questions, please do not hesitate to contact RNA Medical using the contact information below.

RNA Medical, Division of Bionostics, Inc.
 7 Jackson Road
 Devens, MA 01434
 Phone: 978-772-9070 / 800-533-6162
 E-Mail: customerservice@rnamedical.com

BIOMKT-223
 Rev. A