



February 15, 2023

Dear Valued Customer,

CVC223 lot 222341, expiration date 2023-12-31 has an unintentional shift in values on the insert for the level 1 tHb values for the Avox 4000 instrument.

Level 1 lot 25138, expiration date 2024-01-31

Current values

Mean 4.7 g/dL

Range 4.1 – 5.3 g/dL

Corrected values

Mean 5.2 g/dL

Range 4.6 – 5.8 g/dL

The kits in stock will not be reworked to include this update. If you have already purchased this lot, please keep this notification with the kit.

Please note there are no shifts in values for the other levels and no changes to other instruments. If you have any questions, please feel free to reach out using the contact information below.

RNA Medical, Division of Bionostics, Inc.
7 Jackson Road
Devens, MA 01434
Phone: 978-772-9070 / 800-533-6162
E-Mail: customerservice@rnamedical.com

BIOMKT-230 Rev A